



Lp(a)-CONTROL

Lp(a) Control
Latex-Turbidimetry

Quantitative determination of Lipoprotein (a) (Lp(a))

IVD
Store at 2 – 8 °C

PRODUCT CHARACTERISTICS

The control is a lyophilized human serum from human origin containing Lp(a) used to evaluate the precision and the accuracy of Lp(a) determinations in human serum by turbidimetric method.

REAGENTS

Lp(a) CONTROL	Human serum with a low concentration of Lp(a)
---------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence for HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against an Internal Reference Material.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

This control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contamination is prevented during their use. Do not use reagents over the expiration date.

After reconstitution, it is stable for 1 month at 2-8 °C or 3 months at –20 °C.

PROCEDURE

To be used in turidimetric assays.

COMPONENT	VALUE	RANGE	UNITS
Lp(a) CONTROL			mg/dL

Ref: ACC16-028 Cont: 1x1 mL
Lot:
Expiry date:



VITAL DIAGNOSTICS S.r.l. via Balzella 41/G/4, 47100 Forlì, ITALY
Tel: +39 0543 721220 Fax: +39 0543 796001

TF_CTR_LP(a)(01).doc



Lp(a)-CONTROL

Lp(a) Control
Latex-Turbidimetry

Quantitative determination of Lipoprotein (a) (Lp(a))

IVD
Store at 2 – 8 °C

PRODUCT CHARACTERISTICS

The control is a lyophilized human serum from human origin containing Lp(a) used to evaluate the precision and the accuracy of Lp(a) determinations in human serum by turbidimetric method.

REAGENTS

Lp(a) CONTROL	Human serum with a low concentration of Lp(a)
---------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence for HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against an Internal Reference Material.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

This control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contamination is prevented during their use. Do not use reagents over the expiration date.

After reconstitution, it is stable for 1 month at 2-8 °C or 3 months at –20 °C.

PROCEDURE

To be used in turidimetric assays.

COMPONENT	VALUE	RANGE	UNITS
Lp(a) CONTROL			mg/dL

Ref: ACC16-028 Cont: 1x1 mL
Lot:
Expiry date:



VITAL DIAGNOSTICS S.r.l. via Balzella 41/G/4, 47100 Forlì, ITALY
Tel: +39 0543 721220 Fax: +39 0543 796001

TF_CTR_LP(a)(01).doc