



FERR-CONTROL

Ferritin Control

Quantitative determination of FERRITIN

IVD

Store at 2 – 8 °C

PRODUCT CHARACTERISTICS

The control is a synthetic biological fluid containing human ferritin at fixed value diluted in HEPES pH 7.4 buffer containing stabilizers and sodium azide as preservative (<1g/L). It is used to evaluate the precision and the accuracy of ferritin determinations in human serum.

REAGENTS

FERRITIN CONTROL	Synthetic biological fluid containing human ferritin
------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence for HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

PREPARATION

Ready to use

STORAGE AND STABILITY

This control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contamination is prevented during their use. Do not use reagents over the expiration date.

COMPONENT	VALUE	RANGE	UNITS
FERRITIN CONTROL	105	75 -135	ng/mL

Ref: ACC16-030 Cont: 1x2 mL
Lot: 06E04
Expiry date: 11/2007



VITAL DIAGNOSTICS S.r.l. via Balzella 41/G/4, 47100 Forlì, ITALY
Tel: +39 0543 721220 Fax: +39 0543 796001

TF_CTR_FERRITIN(02).doc



FERR-CONTROL

Ferritin Control

Quantitative determination of FERRITIN

IVD

Store at 2 – 8 °C

PRODUCT CHARACTERISTICS

The control is a synthetic biological fluid containing human ferritin at fixed value diluted in HEPES pH 7.4 buffer containing stabilizers and sodium azide as preservative (<1g/L). It is used to evaluate the precision and the accuracy of ferritin determinations in human serum.

REAGENTS

FERRITIN CONTROL	Synthetic biological fluid containing human ferritin
------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence for HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

PREPARATION

Ready to use

STORAGE AND STABILITY

This control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contamination is prevented during their use. Do not use reagents over the expiration date.

COMPONENT	VALUE	RANGE	UNITS
FERRITIN CONTROL	105	75 -135	ng/mL

Ref: ACC16-030 Cont: 1x2 mL
Lot: 06E04
Expiry date: 11/2007



VITAL DIAGNOSTICS S.r.l. via Balzella 41/G/4, 47100 Forlì, ITALY
Tel: +39 0543 721220 Fax: +39 0543 796001

TF_CTR_FERRITIN(02).doc