



μALB-CONTROL

Microalbumin Latex control Turbidimetry

Quantitative determination of microalbumin

IVD
Store at 2 – 8 °C

PRODUCT CHARACTERISTICS

The control is a lyophilized human serum containing human albumin used to evaluate the precision and the accuracy of microalbumin determinations in human urine by turbidimetric method.

REAGENTS

μ-ALB CONTROL	Human serum with a low concentration of microalbumin.
---------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence for HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against the International Reference Material CRM 470/RPPHS (institute for Reference Materials and Measurements, IRMM).

PREPARATION

Reconstitute (→) with 2.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

This control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contamination is prevented during their use. Do not use reagents over the expiration date.

After the reconstitution, it is stable for 10 days at 2-8 °C or 6 months at –20 °C.

PROCEDURE

To be used in turbidimetric assays.

COMPONENT	VALUE	RANGE	UNITS
Microalbumin			mg/L

Ref: ACC16-033 Cont: 1x2 mL
Lot:
Expiry date:



VITAL DIAGNOSTICS S.r.l. via Balzella 41/G/4, 47100 Forlì, ITALY
Tel: +39 0543 721220 Fax: +39 0543 796001

TF_CTR_MICROALBUMIN(01).doc



μALB-CONTROL

Microalbumin Latex control Turbidimetry

Quantitative determination of microalbumin

IVD
Store at 2 – 8 °C

PRODUCT CHARACTERISTICS

The control is a lyophilized human serum containing human albumin used to evaluate the precision and the accuracy of microalbumin determinations in human urine by turbidimetric method.

REAGENTS

μ-ALB CONTROL	Human serum with a low concentration of microalbumin.
---------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence for HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against the International Reference Material CRM 470/RPPHS (institute for Reference Materials and Measurements, IRMM).

PREPARATION

Reconstitute (→) with 2.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

This control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contamination is prevented during their use. Do not use reagents over the expiration date.

After the reconstitution, it is stable for 10 days at 2-8 °C or 6 months at –20 °C.

PROCEDURE

To be used in turbidimetric assays.

COMPONENT	VALUE	RANGE	UNITS
Microalbumin			mg/L

Ref: ACC16-033 Cont: 1x2 mL
Lot:
Expiry date:



VITAL DIAGNOSTICS S.r.l. via Balzella 41/G/4, 47100 Forlì, ITALY
Tel: +39 0543 721220 Fax: +39 0543 796001

TF_CTR_MICROALBUMIN(01).doc