



NEWS

FOR IMMEDIATE RELEASE

Vital Diagnostics Launches the Eon[®] One Semi Automated Photometer with Dilutor

Noble Park, VIC, Australia, September 15, 2010 – Vital Diagnostics' Eon One semi automated photometer has been enhanced with the release of the Eon One dilutor. The use of the Eon One dilutor is convenient, leads to real productivity gains, avoids manual pipetting errors due to incorrect or inaccurate volumes and results in significant improvements in precision when compared to the use of a manual pipette.

The Eon One is a state of the art semi automated photometer that is capable of running all routine clinical chemistry tests. The Eon One dilutor (sold separately) is plugged directly into the Eon one via a dedicated port, and once enabled is ready to use.

The computer embedded into the Eon One instrument controls the dilutor, so aspiration volumes for both reagent and sample for each test parameter programmed into the Eon one are automatically loaded by the system, making it very convenient for the operator.

By simply following easy to read screen prompts, the operator aspirates the required volume of reagent, then sample and dispenses into the reaction tube for the required incubation. The dilutor then prompts the operator to wash before moving to the next sample.

For further information, contact Vital Diagnostics at + 61 2 9894 6988 (ph), + 61 2 9894 6303(fax) or email sales@vitaldiagnostics.com (website www.vitaldiagnostics.com)

Eon and Vital Diagnostics are registered trademarks of Vital Diagnostics.